

# KAZAMA ENT TOTAL CATALOGUE



# KAU-3000

## [HARMONY] - BASIC

### FEATURES

This product is in the scope of medium size of ENT Treatment Unit with new functions to be able to perform treatment and examination in more sanitary environment.


The anti-fog and endoscope holder are designed and mounted on the convenient places to use and the heating device (option) is available to maintain constant temperature of instrument tray.

Unit table is made of sustainable and durable artificial marble and special synthetic resin not only for luxurious looking, but also to prevent any kind of contamination from drug and other chemicals.

### SPECIFICATIONS

|                     |                              |
|---------------------|------------------------------|
| Compressor pump     | 70W, 2.5Kg f/cm <sup>2</sup> |
| Vacuum pump(main)   | Negative pressure 650mmHg    |
| Main suction bottle | 3000cc                       |
| Sub suction bottle  | 1000cc                       |
| Clean system        | Air filter & Regulator       |
| Power consumption   | 600~900W                     |
| Voltage             | 220V                         |
| Dimension           | 880(W) x 620(L) x 740(H)     |
| Weight              | 11.5kg                       |
| Anti-fog device     | 400W                         |

CE 0123



ENT treatment total unit  
KAU-3000 Basic Type

# NET-1500D

## SPECIFICATIONS

- Power Source - AC 220V, 50/60Hz
- Power Consumption - 300VA
- Rotation Angle - 340°
- Return Position - Auto
- Up-Down Stroke - 250(480~730mm)
- Reclining Angle - 90°
- Weight - Approx. 75kg
- Dimension - 600(W) x 730(D) x 1080(H)mm
- Option - Headrest



# NET-1500E

## SPECIFICATIONS

- The Lowest Cost, The Best Cost Effective.
- 45° Reclining Angle.
- Easy to Use
- Noiseless.
- 480mm, The Lowest Height for Seat in the Domestic, Comfortable for Treatment.
- Easy Set Up, can be Handled by One Engineer or Salesperson.
- Exceffect Quality of Practical Application in the Limited Area.





**KAZAMA ENT**

223-612, Seoknam-dong, Seo-gu, Incheon, KOREA, 404-220  
Tel : 82-32-581-1900 Fax : 82-32-582-1309  
<http://www.kazamaent.co.kr> <http://kazamaent.koreasme.com>